The Functional Movement Screen

The system for a simple and quantifiable method of evaluating basic movement abilities

he Functional Movement Screen(FMS) is an innovative system used to evaluate movement pattern quality for clients or athletes. The beauty of the Functional Movement Screen is that a personal trainer, athletic trainer or strength and conditioning coach can learn the system and have a simple and quantifiable method of evaluating basic movement abilities. The FMS only requires the ability to observe basic movement patterns already familiar to the coach or trainer. The key to the Functional Movement Screen is that it consists of a series of simple tests with a simple grading system. The FMS allows a trainer or coach to begin the process of functional movement pattern assessment in individuals without recognized pathology. The FMS is not intended to diagnose orthopedic problems but rather to demonstrate limitations or asymmetries in healthy individuals with respect to basic movement patterns and eventually correlate them with outcomes.

The Functional Movement Screen provides a strength and conditioning coach or personal trainer with an evaluation option that relates closely to what the athlete or client will actually do in training. In a sense, the tests are improved by working on variations of the skills tested. The FMS allows evaluation with tools and movement patterns that readily make sense to both the client and the trainer or coach.

The test is comprised of seven fundamental movement patterns that require a balance of mobility and stability. These fundamental movement patterns are designed to provide observable performance of basic loco motor, manipulative

and stabilizing movements. The tests place the individual in extreme positions where weaknesses and imbalances become noticeable if appropriate stability and mobility is not utilized. It has been observed that many individuals who perform at very high levels during activities are unable to perform these simple movements. These individuals should be considered to be utilizing compensatory movement patterns during their activities, sacrificing efficient movements for inefficient ones in order to perform at high levels. If these compensations continue, then poor movement patterns will be reinforced leading to poor biomechanics.



Gray Cook, MS, PT, OCS, CSCS




Lee Burton, PhD, ATC, CSCS

Test 1: Deep Squat

abduction of the shoulders.

The squat is a movement needed in most athletic events. It is the ready position and is required for most power and lifting movements involving the lower extremities. The deep squat is a test that challenges total body mechanics when performed properly. It is used to assess bilateral, symmetrical and functional mobility of the hips, knees and ankles. The dowel held overhead assesses bilateral, symmetrical mobility of the shoulders as well as the thoracic spine. The ability to perform the deep squat requires appropriate pelvic rhythm, closed-kinetic chain dorsiflexion of the ankles, flexion of the knees and hips and extension of the thoracic spine, as well as flexion and



Test 2: Hurdle Step

The hurdle step is designed to challenge the body's proper stride mechanics during a stepping motion. The movement requires proper coordination and stability between the hips and torso during the stepping motion as well as single leg stance stability. The hurdle step assesses bilateral functional mobility and stability of the hips, knees and ankles. Performing the hurdle step test requires stanceleg stability of the ankle, knee and hip as well as maximal closed-kinetic chain extension of the hip. The hurdle step also requires step-leg open-kinetic chain dorsiflexion of the ankle and flexion of the knee and hip. In addition, the subject must also display adequate balance because the test imposes a need for dynamic stability.

Test 3: In-Line Lunge

This test attempts to place the body in a position that will focus on the stresses as simulated during rotational, decelerating and lateral-type movements. The inline lunge is a test that places the lower extremity in a scissored position, challenging the body's trunk and extremities to resist rotation and maintain proper alignment. This test assesses torso, shoulder, hip and ankle mobility and stability, quadriceps flexibility and knee stability. The ability to perform the in-line lunge test requires stance-leg stability of the ankle, knee and hip as well as apparent closed kineticchain hip abduction. The in-line lunge also requires step-leg mobility of the hip, ankle dorsiflexion and rectus femoris flexibility. The subject must also display adequate stability due to the rotational stress imposed.

Test 6: Trunk Stability Push-up

The actify to disa maintaining straight-leg and gastrating a stable opposite left tive straight hamstring is available. This is difference common required to ity of the open all stability.

The active straight-leg raise tests the ability to disassociate the lower extremity while maintaining stability in the torso. The active straight-leg raise test assesses active hamstring and gastroc-soleus flexibility while maintaining a stable pelvis and active extension of the opposite leg. The ability to perform the active straight-leg raise test requires functional hamstring flexibility, which is the flexibility that is available during training and competition. This is different from passive flexibility, which is more commonly assessed. The subject is also required to demonstrate adequate hip mobility of the opposite leg as well as lower abdominal stability.

Test 5: Active Straight-Leg Raise

The trunk stability push-up tests the ability to stabilize the spine in an anterior and posterior plane during a closed-chain upper body movement. It assesses trunk stability in the sagittal plane while a symmetrical upper-extremity motion is performed. The ability to perform the trunk stability push-up requires symmetric trunk stability in the sagittal plane during a symmetric upper extremity movement. Many functional activities require the trunk stabilizers to transfer force symmetrically from the upper extremities to the lower extremities and vice versa. Movements such as blocking in football and jumping for rebounds in basketball are common examples of this type of energy transfer. If the trunk does not have adequate stability during these activities, kinetic energy will be dispersed, leading to poor functional performance as well as increased potential for micro traumatic injury.





Test 4: Shoulder Mobility

The shoulder mobility screen assesses bilateral shoulder range of motion, combining internal rotation with adduction and external rotation with abduction. It also requires normal scapular mobility and thoracic spine extension. The ability to perform the shoulder mobility test requires shoulder mobility in a combination of motions including abduction/external rotation, flexion/extension and adduction/internal rotation. It also requires scapular and thoracic spine mobility.



Test 7: Rotary Stability

This test is a complex movement requiring proper neuromuscular coordination and energy transfer from one segment of the body to another through the torso. The rotary stability test assesses multi-plane trunk stability during a combined upper and lower extremity motion. The ability to perform the rotary stability test requires asymmetric trunk stability in both sagittal and transverse planes during asymmetric upper and lower extremity movement. Many functional activities require the trunk stabilizers to transfer force asymmetrically from the lower extremities to the upper extremities and vice versa. Running and exploding out of a down stance in football and moving and carrying heavy equipment or objects are examples of this type of energy transfer. If the trunk does not have adequate stability during these activities, kinetic energy will be dispersed, leading to poor performance as well as increased potential for injury.

Scoring the FMS

The individual tests have certain criteria that must be accomplished in order to obtain a high score. The scoring is broken down into four basic criteria: a 3 is given if the individual can perform the movement without any compensations according to the established criteria, a 2 is given if the individual can perform the movement but must utilize poor mechanics and compensatory patterns to accomplish the movement, a 1 is given if the individual cannot perform the movement pattern even with compensations, and finally, a 0 is given if the individual has pain during any part of the movement or test. There are five tests which require bilateral testing; this will result in two scores for those tests. The lowest test score is recorded for the overall score; however, for assessment and data collection purposes, both scores are needed. Three tests: Shoulder Mobility, Trunk Stability Push-up and Rotary Stability have clearing test associated with them that are scored as pass/fail. If a person fails this part of the test, then a 0 is given as the overall score.

The FMS is an assessment technique, which attempts to identify imbalances in mobility and stability during fundamental movement patterns. This assessment tool is thought to exacerbate the individual's compensatory movement problems, allowing for easy identification. It is these movement flaws that may lead to breakdown in the kinetic linking system, causing inefficiency and microtrauma during activity.

The FMS should be introduced as part of the preplacement/pre-participation physical exam-ination to determine deficits that may be overlooked during the traditional medical and performance evaluations. In many cases, muscle flexibility and strength imbalances along with previous injuries may not be identified. These problems, which have been acknowledged as significant risk factors for injury, will be identified using the FMS. This movement-based assessment will pinpoint functional deficits related to proprioceptive, mobility and stability weaknesses. If these risk factors can be identified and addressed utilizing the FMS, then decreases in injuries and improved performance should follow.